

# Analytical characteristics of new Faraday amplifier boards equipped with 1e<sup>11</sup> and 1T resistors for use in high precision TIMS isotope ratio measurements.

Excellence in mass spectrometry

# Poster number V23B-2808



#### Introduction

Thermal Ionization Mass Spectrometers (TIMS) are capable of generating the highest precision and accuracy in isotope ratio determinations. Low mass bias, accurate mass fractionation corrections, high sensitivity and the absence of isobaric interferences are key contributions to a TIMS performance. However, the ability to quantitatively measure ion currents in the picoamp range with minimal calibration are essential ingredients to a precise isotope ratio measurement.

We present isotope ratio data obtained from the Phoenix TIMS equipped with new **X**act amplifier boards, which demonstrates the current performance levels of a TIMS.

#### Xact amplifier boards.

The new Xact amplifier boards are mounted within the evacuated and peltier cooled amplifier housing on the Phoenix TIMS. The boards are also back compatible with IsoProbe-T, IsoProbe-P and Sector 54 mass spectrometers.

The board material itself has been designed to be hydrophobic as possible as the Figure 1. Gain measurements on 9 boards over 30 days presence of residual water in the board material which degasses over time is probably a principal cause of gain instability over time. The resistors have a 25ppm/°C **Resistor/Amplifier response (speed)** temperature coefficient which is almost 10 times lower than the previous type of re-Figure 2 shows the response of a 1e<sup>11</sup> ohm resistor to removal of an 8e<sup>-11</sup> amp ion sigsistors which have been in use since the early 1990's. Since temperature regulation nal. The Y co-ordinate shows the intensity in ppm of signal, and the x axis, time in within the amplifer bin is kept to 16°C+/-0.02°C then the theoretical change in gain seconds. After 0.5 seconds the intensity is less than 10ppm of the signal, i.e, 8e<sup>-16</sup>A. would only be 0.5ppm. The amplifier bin is permanently pumped using an oil free After less than 1 second it is at baseline. scroll pump to minimize any ingress of oil from the pump, and maintain a vacuum of <0.1mbar.

### **Faraday Noise**

The Noise levels for different integration times of the new boards are shown in Table 1, for both standard 1e<sup>11</sup> ohm resistor boards and for boards equipped with 1T resistors. In theory the 1T boards should be 3 times quieter than the 1e<sup>11</sup> ohm boards, but this is not the case. The 1T boards are typically 1.4 to 1.8 times quieter than the 1e<sup>11</sup> ohm boards, with the biggest benefit occurring for 10 second integrations.

Looking at it another way the noise characteristics for the 1T resistor for a 10 second integration is the same as a 1e<sup>11</sup>ohm resistor for a 40 second integration. In effect for small ion beams you would have to integrate the baseline or the signal four times longer for the 1e<sup>11</sup> ohm resistors than for the 1T resistors. This is clearly where the 1T resistors are advantageous.

Zenon Palacz (zenon.palacz@isotopx.com), Isotopx Ltd., Millbrook Court, Midpoint 18, Middlewich, Cheshire CW10 0GE, United Kingdom

Noise (Amps)						
Integration time (seconds)	1e <sup>12</sup> ohm n=5 (amps)	1e <sup>11</sup> ohm n=4 (amps)	<b>1e<sup>11</sup>/1e<sup>12</sup></b>			
1	2.3E-16	3.8E-16	1.7			
5	1.0E-16	1.7E-16	1.6			
10	6.5E-17	1.2E-16	1.8			
20	5.1E-17	7.3E-17	1.4			
40	3.6E-17	6.3E-17	1.7			
60	3.6E-17	5.1E-17	1.4			
300	3.5E-17	4.9E-17	1.4			

**Table 1.** Noise comparison between 1T and  $1e^{11} \Omega$  resistors.

# Gain stability

Figure 1 shows gain measurements measured on 9 boards (8 relative to the axial board).

The measurements were made over 25 days on an instrument in a non temperature controlled environment which experienced a minimum of 10°C diurnal temperature variation.

The gains (with 1 exception, series2) vary +/-5ppm over this period. The gains appear to vary in a systematic manner, which is highlighted by the excursion at day 12. All the boards seem to show the same variation at the same time. This is variation is probably due to external temperature variations in the environment, and quite possibly due to variations in the stability of the reference voltage which is supplied by the FAC unit which is not temperature regulated.





*Figure 2.* Response of 1e<sup>11</sup> ohm resistor after removal of an 8e<sup>-11</sup>A beam

For applications which require fast peak jumping e.g. high precision multidynamic measurements, there is virtually no residual signal within 1 second of a mass jump. For comparison the previous version boards would attain 5ppm of signal within 2 seconds, and the pre-1990 version produced 10ppm of signal within 2 seconds.

Figure 3 shows the response of 1T resistors on the same amplifier board. The figure shows the trace for several resistors, the response is the same. Following an initial overshoot in the first second, there is a rebound, before finally returning to baseline



*Figure 3.* Response of 1T ohm resistor after removal of an 8e<sup>-11</sup>A beam

(<5ppm of signal) in 4 seconds. This is approximately a factor of ten faster than the previous design of boards, and makes the possibility of multidynamic measurements possible with these boards.

## Sr isotope ratio data

High Precision multidynamic measurements would be expected from low noise, fast resistors. Figure 4 is an example of Sr using a 5 volt <sup>88</sup>Sr ion beam. The stability of the *Figure 7. Replicate measurements of 2ng NdO*<sup>+</sup> gain between the detectors is demonstrated by static Sr isotope ratio measurements shown in Figure 5. There is no drift in the data which would indicate a drift in the gain Even for small samples of Sr (<1ng) 1e<sup>11</sup> ohm resistors are adequate in producing high (or a collector efficiency change). The oscillatory nature of the data is probably due to precision isotope ratios (Table 2), simply because the ionization efficiency of Sr is very how many baseline measurements were made during the analysis. high with TaF activator (>15% for sub nanogram levels)



*Figure 4.*<sup>87</sup>Sr/<sup>86</sup>Sr data acquired using Xact 1e<sup>11</sup> ohm boards.



Figure 5. Gain stability between detectors

# Small sample measurements

The ability to measure small ion signals is shown in Figure 6. Here a Nd ion beam is reduced in size. The errors expand as the ion beam is reduced as would be expected, however, both 1T and 1e<sup>11</sup> ohm boards remain coherent down to 1e<sup>-13</sup>A or 10mv of

<sup>142</sup>Nd or 4mv of <sup>143</sup>Nd. The errors expand below this level.



*Figure 6. Precision as Nd ion beam intensity is reduced* 

For small sample analyses of Nd, analysis using NdO<sup>+</sup> is preferable as the ionization efficiency is about 20% which is ten times higher than Nd<sup>+</sup>. Replicate measurements of NdO<sup>+</sup> show that for 2 ng of Nd, there is very little difference between 1e<sup>11</sup> ohm and 1T resistors, simply because the ion intensities are too high ~2-4e<sup>-12</sup>A, which is consistent with Figure 7.



	Total Evaporation NBS 987				
	100pg	200pg	500pg	1ng	
	0.710244	0.710320	0.710249	0.710266	
	0.710255	0.710279	0.710263	0.710262	
	0.710285	0.710269	0.710265	0.710274	
	0.710242	0.710222	0.710264	0.710278	
	0.710256	0.710298	0.710273	0.710253	
	0.710320	0.710229	0.710253	0.710265	
	0.710299	0.710250	0.710249	0.710257	
	0.710229		0.710249	0.710273	
	0.710285		0.710272	0.710254	
	0.710303		0.710268	0.710258	
mean	0.710272	0.710267	0.710261	0.710264	
1SD	0.000031	0.000036	0.000010	0.000009	
1RSD (PPM)	43	50	14	12	

Table 2. Sr isotope ratio measurements with 1e<sup>11</sup>ohm resistors.

### Conclusions

- New **X**act amplifier boards are quiet and fast.
- Gain calibrations are stable enough not to require gain calibration be tween samples or the use of dynamic amplifiers.
- Noise levels of 1T resistors are less than 2 times better than 1e<sup>11</sup> ohm resistors and do not produce significant benefits for nanogram levels of Nd or >100pg of Sr due to the high ionization efficiency of these elements.
- 1T resistors can probably prove advantageous for very small ion beams (<1e<sup>-14</sup>A) and for transient ion signals.